

Silicone Instrument Mats

National Surgical Corporation Silicone Instrument Mats: high quality, large range, protect instruments

HOLDS INSTRUMENTS FOR STERILISATION

Place delicate and very delicate instruments on mats to ensure no movement during the sterilisation and cleaning process

DELICATE OR VERY DELICATE INSTRUMENTS

Different sizes available to hold delicate instruments with 8mm distance between pins and very delicate instruments with 3mm between pins or 19mm pins for enhanced airing of instruments

DIFFERENT SIZES AVAILABLE

Large range of sizes available and most can be cut to suit tray/basket to securely hold various instruments

CAN FIT IN STERILISATION BASKETS

Ideal for inserting into any type of sterilisation tray or basket, Medicon Medicontainers. Mats can be cut to size to suit tray

PROTECTS INSTRUMENTS

Silicone mats hold instruments securely to ensure they won't move, fall or accidentally be lost during cleaning and sterilisation

STERILISATION BASKETS AVAILABLE

Large range of sterilisation baskets, trays, micro mesh baskets and Medicontainers also available




National Surgical Corporation
is certified to
AS/NZS ISO 9001:2015



BE SUPPORTED BY AN
AUSTRALIAN COMPANY
PROVIDING VALUE BEYOND THE PRODUCT

Silicone Instrument Mats



Product Code	Description
For delicate instruments	
49.97.60	Silicone mats for the fixation of delicate instruments - 6mm, 460x240mm
49.97.61	Silicone mats for the fixation of delicate instruments - 6mm, 275x125mm
49.97.62	Silicone mats for the fixation of delicate instruments - 6mm, 240x220mm
49.97.63	Silicone mats for the fixation of delicate instruments - 6mm, 260x250mm
49.97.64	Silicone mats for the fixation of delicate instruments - 6mm, 265x115mm
49.97.65	Silicone mats for the fixation of delicate instruments - 13mm, 460x240mm
49.97.66	Silicone mats for the fixation of delicate instruments - 13mm, 275x125mm
49.97.67	Silicone mats for the fixation of delicate instruments - 13mm, 240x220mm
49.97.68	Silicone mats for the fixation of delicate instruments - 13mm, 260x250mm
49.97.75	Silicone mats for the fixation of delicate instruments - 8mm, 460x240mm
49.97.76	Silicone mats for the fixation of delicate instruments - 8mm, 275x125mm
49.97.77	Silicone mats for the fixation of delicate instruments - 8mm, 240x220mm
49.97.78	Silicone mats for the fixation of delicate instruments - 8mm, 260x250mm
49.97.80	Silicone mats for the fixation of delicate instruments - 8mm, 467x240mm
49.97.81	Silicone mats for the fixation of delicate instruments - 8mm, 265x115mm
49.97.82	Silicone mats for the fixation of delicate instruments - 13mm, 467x240mm
49.97.83	Silicone mats for the fixation of delicate instruments - 13mm, 265x115mm
49.97.87	Silicone mats for the fixation of delicate instruments - 6mm, 545x262mm
49.97.88	Silicone mats for the fixation of delicate instruments - 6mm, 490x240mm
For very delicate instruments	
49.97.70	Silicone mats for the fixation of very delicate instruments - 3mm, 460x240mm
49.97.71	Silicone mats for the fixation of very delicate instruments - 3mm, 275x125mm
49.97.72	Silicone mats for the fixation of very delicate instruments - 3mm, 240x220mm
49.97.73	Silicone mats for the fixation of very delicate instruments - 3mm, 260x250mm
49.97.74	Silicone mats for the fixation of very delicate instruments - 3mm, 265x115mm
49.97.85	Silicone mats for the fixation of very delicate instruments - 3mm, 545x262mm
49.97.86	Silicone mats for the fixation of very delicate instruments - 3mm, 490x240mm
Multipurpose silicone mats	
SH18444	Silicone pin mat - 33cm x 46cm
SH18445	Silicone pin mat - 33cm x 46cm
ACC-5015	Silicone mat - 515 x 240 x 15mm
ACC-5016	Silicone mat - 545 x 565 x 19mm

National Surgical Corporation has a full range of surgical instruments and equipment please visit our Surgical Instrument website: www.surgicalinstrument.net.au

CUSTOMER SERVICE

Phone: 1800 138 138 (free call)

Fax: 1800 138 140 (free fax)

info@nationalsurgical.com.au

www.nationalsurgical.com.au

Contact your National Surgical Corporation Account Manager for more information



NATIONAL SURGICAL CORPORATION

A proudly founded and owned Australian company - providing value to our customers beyond the product