



# Vacuform 2.0 Surgical Bedding

## *National Surgical Corporation Vacuform Surgical Bedding: the feeling of weightlessness in surgery*

### **INNOVATIVE MICROSPACE SYSTEM**

Microspace system allows the mattress to adapt perfectly to the body using microballs of styrofoam, stabilised by vacuum to give a lasting, completely supportive form

### **REDUCE RISK OF BED SORES**

Reduces the risk of bed sores (Decubitus prophylaxis) during lengthy procedures or in the elderly by adapting to the patient's body resulting in an even distribution of contact pressure, preserving micro-circulation

### **SAFETY DURING SURGERY**

The patient can be fixed quickly and safely in any required position, giving the surgeon optimum access to the operation site

### **LARGE RANGE OF MATTRESSES**

Large range of mattresses for operations to the shoulders, hips, legs, feet and special designs for lateral positions. Any size and shape can be developed to meet your requirements

### **MADE FROM SLIP PROOF, LATEX FREE PVC**

Mattresses are made from special, soft PVC sheeting, making them slip proof and stable

### **PLASTIC VALVE**

Allows for one handed operation

### **SURGICAL HOOD - PERMEABLE TO X-RAYS**

Specially developed hood envelops the whole head and can be used to immobilise the patient's head during CT scans as it's permeable to X-rays

### **SURGICAL CUSHIONS**

Range of surgical cushions also adapt to the patient's head for superior comfort. Full X-ray, MRI and CT compatibility

### **ACCESSORIES AVAILABLE**

Range of surgical bedding accessories including surgical sheets, non-slip sheeting and patient positioning aids

# Vacuform 2.0 Surgical Bedding



MicroSpace filling in the evacuated state      Conventional filling in the evacuated state

Image	Product Code	Description
<b>Vacuform Surgical Mattress</b>		
	3466105	Vacuform Universal surgical mattress Size: 75x200cm
	3466205	Vacuform Universal surgical mattress Size: 65x200cm
	3466104	Vacuform Universal surgical mattress Size: 100x200cm
	3466505	Vacuform Lower Leg & Foot surgical mattress - Size 75x200cm
	3466605	Vacuform Shoulder surgical mattress with separate cushion Size 75x200cm
	3466605-k	Vacuform Shoulder surgical mattress with separate cushion (30x35cm) Size 75x200cm
	3466305	Vacuform Side position surgical mattress according to Murpf Size 95/70 x 212cm
	3466405	Vacuform Shoulder surgical mattress with Head Positioning Size 75x200cm
	3465205	Vacuform Short surgical mattress Size: 75x125cm
	3466005	Vacuform Short surgical mattress Size: 100x100cm
	3465005	Vacuform Short surgical mattress Size: 75x100 cm
	3465105	Vacuform Short surgical mattress with cutout - Size: 75x100cm

Image	Product Code	Description
	3465305	Vacuform Short surgical mattress with cutout - Size: 75x125cm
	3466006	Vacuform Short surgical mattress with cutout - Size: 100x100cm
	3466007	Vacuform Side position Short surgical mattress during a minimally invasive hip operation - Size: 95/80x105cm
	3466008	Vacuform Plasterbed, reusable Size: 75x200cm
<b>Vacuum cushions</b>		
	3464005	Vacuum Cushion - Positioning of back, infants - Size 60x80cm
	3463005	Vacuum Cushion - Positioning of leg Size 40x70cm
	3462005	Vacuum Cushion - Positioning of arm Size 30x70cm
	3461005	Vacuum Cushion - Positioning of head and hand - Size 30x40cm
	3467005	Vacuum Cushion - Neck Positioning Size 35x40cm
<b>Vacuform Surgical Hood</b>		
	3464700	Vacuform surgical hood - small
	3464701	Vacuform surgical hood - large

Range of accessories also available including cushions made from rubber and synthetic for lumbar, shoulder, head, pelvis and ulnar

**CUSTOMER SERVICE**  
 Phone: 1800 138 138 (free call)  
 Fax: 1800 138 140 (free fax)  
 info@nationalsurgical.com.au  
 www.nationalsurgical.com.au

Contact your National Surgical Corporation Account Manager for more information



**NATIONAL SURGICAL CORPORATION**

A proudly founded and owned Australian company - providing value to our customers beyond the product